


Wireless Portable X-Ray **Biox**

Anytime
Anyplace



Handy and Powerful

Excellent Image

Ergonomic Design



www.digimed.co.kr



Total weight 2Kg, 3mA X-ray Tube current

- Simple operation.
- Wide graphic LCD.
- Light weight and ergonomic design.
- Highest image quality and lowest patient dosage.
- Zero leakage radiation
- High frequency, constant-potential X-ray generator.
- Optimized to analog film and digital sensor both
- Easy work with optional arms.



BIOX Specification

X-Ray Generator	DC High Frequency	Biox
X-Ray Insulation	Tube Voltage	60kV
	Tube Current	3mA
X-Ray Tube	Anode Type	Stationary
	Anode Angle	20°
	Focal spot size	0.8mm
	Heat capacity	8.5Khu
Exposure Time	0.01~1.6 sec(0.01 step)	
Total Filtration	1.6mm Al	
SSD	100mm, 210mm(Long cone)	
Total Weight	2.2kg	
Dimension	180(W) x 73(D) x 138(H) mm	
Power Requirement	Battery	DC 24V
	Battery Charger	AC 100~240V
	Power Supply(Stand)	AC 100~240V



Manufacture



#311, ACE HIGH-END TOWER 3, 371-50,
Gasam-Dong, Geumcheon-Gu, Seoul, Korea
Tel.: 82-2-2624-1551 Fax.: 82-2-2624-1553
<http://www.digimed.co.kr>
E-mail : digi-med@hanmail.net

Distributed by

