





Portable X-ray System

MÍNÍX

Anytime..

Anyplace..

Smaller and Lighter Design
High Frequency X-ray
Low Radiation Dosage



DigiMed

www.digimed.co.kr

New Generation Portable X-ray System

- More compact and smart
- Wide graphic LCD control
- High Frequency X-ray
- Low patient/operator dosage
- Simple operation



Specification

Model Name		MINIX-S	MINIX-V
X-Ray Installation	Tube Voltage	60 kV	70 kV
	Tube Current	2 mA	2 mA
X-Ray Tube	Anode Type	Stationary	
	Anode Angle	20°	12.5°
	Focal Spot	0.8 mm	0.4 mm
	Heat Capacity	8.5 KHU	1400J
Exposure Time Setting		0.01 ~ 2.00 Sec. (0.01 Sec. step)	
Total Filtration		1.8 mmAl.	
SSD		100/200 mm	200 mm
Total Weight		1.7 kg	1.9 kg
Dimension		123(w) x 130(d) x 134(h)	123(w) x 220(d) x 134(h)
Power Requirement	Battery	DC 24V	
	Charger	AC 100 ~ 240 V	

Manufactured by



DigiMed

Distributed by

309~311, 318-ho, 145

Gasam digital 1-ro, Geumcheon-gu, Seoul, Korea

Tel.: +82-2-2624-1551 Fax.: +82-2-2624-1553

digi-med@hanmail.net

www.digimed.co.kr