

Portable X-ray System

PROX



Anytime!
Anyplace!

- ◆ User-friendly design
- ◆ High frequency X-ray
- ◆ Fine image but Low dose radiation
- ◆ Anatomical select buttons for Easy Setting
- ◆ Working great with Film and Sensor



DigiMed

Portable X-ray System, **PROX**



- Light weight & ergonomic design
- Simple operation
- Excellent image quality
- Lowest patient & operator dosage
- Wide graphic LCD
- High frequency, constant-potential x-ray generator
- Easy control with anatomical select buttons



Specification

Model Name		PROX
X-Ray Generator		DC High Frequency
X-Ray Installation	Tube Voltage	60 kV
	Tube Current	2 mA
X-Ray Tube	Anode Type	Stationary
	Anode Angle	20°
	Focal Spot	0.8 mm
	Heat Capacity	8.5 KHU
Exposure Time Setting		0.01 ~ 1.60 Sec. (0.01 Sec step)
Total Filtration		1.8 mmAl
SSD		100 mm / 200 mm(Long cone)
Total Weight		2.2 kg
Dimension		180(w) x 145(d) x 138(h)
Power Requirement	Battery	DC 24 V
	Charger	AC 100 ~ 240 V

Manufactured by



DigiMed

309~311, 318-ho, 145

Gasam digital 1-ro, Geumcheon-gu, Seoul, Korea

Tel.: +82-2-2624-1551 Fax.: +82-2-2624-1553

www.digimed.co.kr

E-mail : digi-med@hanmail.net

Distributed by

