


# Wireless Portable X-Ray **DÍOX**



**Compact size and light weight**

**Excellent image quality**

**Wireless remote controller for X-ray exposure**

**DIGIMED**

[www.digimed.co.kr](http://www.digimed.co.kr)



## Total weight 1.8Kg, 2mA X-ray tube current

- Wide graphic LCD.
- Highest quality image and lowest patient dosage.
- High frequency operated constant potential X-ray generator.
- Pre-programmed quick setting makes the operation fast and easy.
- Safe and powerful Li-Polymer rechargeable battery with PCM.
- Overlapped internal radiation shield is specially designed to protect the operator from radiation.
- Mobile stand, wall mount stand.



## Specification

X-Ray Generator	DC High Frequency	DIOX-602
X-Ray Installation	Tube Voltage	60kV
	Tube Current	2mA
X-Ray Tube	Anode Type	Stationary
	Anode Angle	20°
	Focal Spot	0.8mm
	Heat Capacity	8.5KHU
Exposure time Set	0.01~1.60sec(0.01Sec step)	
Total Filtration	1.6mm Al	
SSD	200mm(Long cone)	
Total Weight	1.8Kg	
Dimension	135(W) x 150(D) x 175(H)mm	
Power requirement	Battery	DC 24V
	Battery Charger	AC100~240V

Manufacture

**DIGIMED**

401, 402, 406~411-ho,  
101, Gasan digital 2-ro, Geumcheon-gu  
Seoul 08505 Rep. of Korea  
T.82-2-6265-1551 F.82-2-6265-1553  
www.digimed.co.kr  
digi-med@hanmail.net

Distributed by

