

Portable X-ray System

# MiniX




*Anytime..*

*Anyplace..*

Smaller and Lighter Design

High Frequency X-ray

Low Radiation Dosage



**DIGIMED**

[www.digimed.co.kr](http://www.digimed.co.kr)



# New Generation Portable X-ray System

- More compact and smart
- Wide graphic LCD control
- High Frequency X-ray
- Low patient/operator dosage
- Simple operation



## Specification

| Model Name            |               | MINIX-S                           | MINIX-V                  |
|-----------------------|---------------|-----------------------------------|--------------------------|
| X-Ray Installation    | Tube Voltage  | 60 kV                             | 70 kV                    |
|                       | Tube Current  | 2 mA                              | 2 mA                     |
| X-Ray Tube            | Anode Type    | Stationary                        |                          |
|                       | Anode Angle   | 20°                               | 12.5°                    |
|                       | Focal Spot    | 0.8 mm                            | 0.4 mm                   |
|                       | Heat Capacity | 8.5 KHU                           | 1400J                    |
| Exposure Time Setting |               | 0.01 ~ 2.00 Sec. (0.01 Sec. step) |                          |
| Total Filtration      |               | 1.8 mmAl.                         |                          |
| SSD                   |               | 100/200 mm                        | 200 mm                   |
| Total Weight          |               | 1.7 kg                            | 1.9 kg                   |
| Dimension             |               | 123(w) x 130(d) x 134(h)          | 123(w) x 220(d) x 134(h) |
| Power Requirement     | Battery       | DC 24V                            |                          |
|                       | Charger       | AC 100 ~ 240 V                    |                          |

Manufactured by

**DIGIMED**

401, 402, 406~411-ho,  
101, Gasan digital 2-ro, Geumcheon-gu  
Seoul 08505 Rep. of Korea  
T.82-2-6265-1551 F.82-2-6265-1553  
www.digimed.co.kr  
digi-med@hanmail.net

Distributed by

