

PROX-S

Portable X-ray System



- ◆ TFT LCD *Touch Screen*
- ◆ *High frequency X-ray*
- ◆ Fine image with *Low dose* radiation
- ◆ Anatomical select buttons for *Easy Setting*
- ◆ Working great with *Film and Sensor*

DIGIMED

www.digimed.co.kr

Portable X-ray System, **PROX-S**



- Wide Touch Screen
- Internal Radiation Shield
- User-friendly Design
- Low Patient & Operator Dosage
- Easy Control & Simple Operation



- *Optional Mounting Arms
- Mobile stand arm
 - Wall mount arm
 - Portable stand arm

Specification

Model Name		PROX-S
X-Ray Generator		DC High Frequency
X-Ray Installation	Tube Voltage	60 kV
	Tube Current	2 mA
X-Ray Tube	Anode Type	Stationary
	Anode Angle	20°
	Focal Spot	0.8 mm
	Heat Capacity	8.5 KHU
Exposure Time Setting		0.01 ~ 2.00 Sec. (0.01 Sec. step)
Total Filtration		2.8 mm Al.
SSD		100 mm / 200 mm (Long cone)
Total Weight		2.2 kg
Dimension		180(w) x 145(d) x 138(h)
Power Requirement	Battery	DC 24 V
	Charger	AC 100 ~ 240 V

Manufactured by

DIGIMED

401, 402, 406~411-ho,
101, Gasan digital 2-ro, Geumcheon-gu
Seoul 08505 Rep. of Korea
T.82-2-6265-1551 F.82-2-6265-1553
www.digimed.co.kr
E-mail : digi-med@hanmail.net

Distributed by

