

Portable X-ray System
DHX-70H

Safe & Easy!



DIGIMED

Technical Specification

DIGIMED

- Tube Voltage : 70 Kv
- Tube Current : 2 mA
- Focal Spot : 0.4 mm
- Time Range : 0.05 - 1 sec
- Total Filtration : Over 1.5 mmAl
- SSD : 200 mm
- Weight : 1.9 kg
- Dimension : 168 x 228 x 118 mm
- Power Source : 14.8V Rechargeable Battery



**Switchable
Rectangular
Emission cone
&
Backscatter
Shield**



RELIABLE X-RAY SOURCE
Meet Safety Requirements



EASE OF USE
Optimal Grip in your hand
Patient-friendly



ROTARY SWITCH
Easy and Quick Setup




SAFETY PROTECTION
Internal Radiation Protection



WIDE TFT COLOR LCD
Clearer readability

*Just right where
you need it!*



Manufactured by,

DIGIMED Co., Ltd.

401, 402, 406~411-ho, A-dong
101, Gasan digital 2-ro, Geumcheon-gu
Seoul 08506 Rep. of KOREA
TEL. 82.2.2624.1551 FAX. 82.2.2624.1553
digimed@digimed.co.kr

www.digimed.co.kr

Distributed by,

