

Portable X-ray System  
**HYBRID Pro**

*Safe & Easy!*



**DIGIMED**



# Technical Specification

# DIGIMED

- Tube Voltage : 70 Kv
- Tube Current : 2 mA
- Focal Spot : 0.4 mm
- Time Range : 0.05 - 1 sec
- Total Filtration : Over 1.5 mmAl
- SSD : 200 mm
- Weight : 1.9 kg
- Dimension : 168 x 228 x 109 mm
- Power Source : 14.8V Rechargeable Battery



**Switchable  
Rectangular  
Emission cone  
&  
Backscatter  
Shield**



**RELIABLE X-RAY SOURCE**  
Meet Safety Requirements



**EASE OF USE**  
Optimal Grip in your hand  
Patient-friendly




**ROTARY SWITCH**  
Easy and Quick Setup



**SAFETY PROTECTION**  
Internal Radiation Protection



**WIDE TFT COLOR LCD**  
Clearer readability



*Just right where  
you need it!*

Manufactured by,

**DIGIMED Co., Ltd.**

401, 402, 406~411-ho, A-dong  
101, Gasan digital 2-ro, Geumcheon-gu  
Seoul 08505 Rep. of KOREA  
TEL. 82.2.6265.1551 FAX. 82.2.6265.1553  
digimed@digimed.co.kr  
[www.digimed.co.kr](http://www.digimed.co.kr)

Distributed by,