

HYBRID

PORTABLE X-RAY SYSTEM

***Remember.
Bring our care for you; Safety and Quality.***

HYBRID S70



DIGIMED

Technical Specification

- Tube Voltage : 70 Kv
- Tube Current : 2 mA
- Focal Spot : 0.4 mm
- Time Range : 0.05 - 1 sec
- Total Filtration : Over 1.5 mmAl
- SSD : 200 mm
- Weight : 1.9 kg
- Dimension : 146 x 227 x 126 mm
- Power Source : 14.8V Rechargeable Battery

DIGIMED



*Switchable
Rectangular
Emission cone*



RELIABLE X-RAY SOURCE
Safety Requirements cleared



EASY TO HANDLE
Best Grip on your palm
Harmony with Patient



ROTARY SWITCH
Easy and Fast Setting




STRONGER SAFEGUARD
Internal Radiation Protection



WIDE TFT COLOR LCD
More clearly visible

*Just right there
you need it!*



Manufactured by,
DIGIMED Co., Ltd.

401, 402, 406~411-ho, A-dong,
101, Gasan digital 1-ro, Geumcheon-gu
Seoul 08505 Rep. of KOREA
TEL. 82.2.6265.1551 FAX. 82.2.6265.1553
digimed@digimed.co.kr
www.digimed.co.kr

Distributed by,